

NEW STANDARD OF RESIDUAL WHITE BLOOD CELL COUNTER

# Walk-away!

ADAM rWBC HT

50 min / 50 tests



High Throughput  
residual White Blood Cell Counter

# ADAM rWBC HT

Walk-away residual White Blood Cell Counter

ADAM-rWBC HT gives you walkaway convenience and significantly improves blood bank's productivity. It takes 50 minutes to count 50 samples.

ADAM-rWBC HT counts residual white blood cells (rWBCs) with an unparalleled ease and speed.



**Not required** ADAM-rWBC HT, eliminates the need for:

Warm-up

Daily maintenance


Professional operator

Vortex

Dark room incubation

Gating adjustment

## Features




### **FAST** measurement

ADAM-rWBC HT counts up to 50 samples at a time, taking only **45 seconds** per sample.



### **EASY** to use

ADAM-rWBC HT process is simple and easy for multi counting.





### **ACCURATE** and reliable result

- $r^2=0.989$
- Substantially equivalent to flow cytometry

# Principle


ADAM-rWBC HT stains nuclei of White Blood Cell using Propidium Iodide (PI) and detects its fluorescent expression. It automatically focuses on the slide and counts the cells from 63 images, then averages out the counting results to increase the accuracy and reliability.




# Performance | Comparison to flow cytometry

Refer to the following comparison of residual blood cell enumeration methods between flow cytometry (Leucocount) and ADAM-rWBC series in using different amounts of white blood cells aliquots. This comparison was performed using both RBC and platelet samples.

## RBC product



## Platelet products



## Precision test

Stain-to-Stain Precision - RBCs							
WBC / $\mu$ L target	Site#	Mean	SD	Total CV%	Mean	SD	Total CV%
		Unit 01			Unit 02		
0-1	1	<1	0.43	NA	1	0.47	44.15
	2	<1	0.62	NA	1	0.96	69.91
	3	<1	0.19	NA	<1	0.60	NA
5-10	1	6	1.05	17.20	8	0.89	10.84
	2	6	1.85	28.97	10	1.32	13.49
	3	7	0.62	9.39	6	0.98	17.57
20-30	1	26	2.11	7.98	27	2.34	8.53
	2	20	0.89	4.38	24	2.33	9.70
	3	25	2.15	8.49	27	3.60	13.59
50-60	1	49	2.77	5.64	54	3.26	6.07
	2	47	3.47	7.42	64	4.59	7.21
	3	54	2.39	4.41	55	3.67	6.71
80-100	1	82	1.91	2.33	90	7.37	8.19
	2	73	6.35	8.69	84	5.14	6.09
	3	91	4.27	4.68	89	4.66	5.24

Stain-to-Stain Precision - Platelets							
WBC / $\mu$ L target	Site#	Mean	SD	Total CV%	Mean	SD	Total CV%
		Unit 01			Unit 02		
0-1	1	1	0.81	60.05	<1	NA	NA
	2	<1	0.22	NA	<1	0.42	NA
	3	<1	0.31	NA	<1	0.44	NA
5-10	1	8	0.85	10.77	8	0.83	11.02
	2	5	1.08	20.66	9	0.87	9.54
	3	6	0.92	14.84	6	1.17	18.37
20-30	1	26	1.79	7.00	25	2.49	10.07
	2	15	2.08	13.89	26	1.64	6.34
	3	28	2.65	9.47	29	2.02	6.88
50-60	1	52	2.84	5.48	55	2.74	4.96
	2	32	2.57	7.94	49	1.68	3.42
	3	61	3.07	5.04	61	3.99	6.58
80-100	1	92	5.89	6.39	90	3.93	4.38
	2	53	2.27	4.30	84	3.03	3.61
	3	98	3.39	3.45	100	3.69	3.67

\* American Red Cross Holland Lab, Rockville, MD

## ADAM-rWBC HT

Cat No. ADAM-rWBC HT

Analysis time	45 sec/test (50 min/50 tests)
Measuring range	1~100 cells / $\mu$ L
Voltage	100-240V~
Current consumption	5.0A
Weight	120 Kg
Size (WxDxH)	1020 x 670 x 730 mm



## ADAM-rWBC Kit

Cat No. ADHK-050

Standard bead	7 mL
r-Solution	25 mL
r-Slide	50 pcs (Loading sample volume : 100 $\mu$ L/test)



### Not affected by temperature

Standard bead and r-Solution can be used without being affected by temperature variations and should not be negatively impacted by shipping.



## Ordering information

Catalog number	Product name
ADAM-rWBC HT	ADAM-rWBC HT

Catalog number	Product name
ADHK-050	ADAM-rWBC Kit (50 tests)



Scan for more information.



### NanoEntek, Inc.

851-14, Seohae-ro, Paltan-myeon, Hwaseong-si, Gyeonggi-do, 18531, Korea  
Tel : +82-2-6220-7940 / Fax : +82-2-6220-7999

### NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA  
Tel : +1-781-472-2558 / Fax : +1-781-790-5649

website [www.nanoentek.com](http://www.nanoentek.com)

e-mail [ivdst@nanoentek.com](mailto:ivdst@nanoentek.com)