
FREND PCT

REF FRPC 025

- English
- German
- Italian
- Spanish
- French
- Portuguese

NESPI-FRPC-001ML(V.0.1)

FREND PCT

Quantitative assay for Procalcitonin

FRPC 025

For *in vitro* diagnostic use only

1. Intended use

The FREND PCT test is designed for *in vitro* diagnostic use only for the quantitative measurement of Procalcitonin (PCT) in human serum. The FREND PCT microfluidic flow cartridge is designed for use in the FREND™ System fluorescent immunoassay reader. For professional use only.

2. Summary and explanation of test

PCT is a precursor of the hormone calcitonin and is produced physiologically by thyroid C cells. PCT concentrations are usually very low in healthy people but is secreted in various extrathyroidal neuroendocrine tissues by circumstances of systemic inflammation, particularly bacterial infection. Concentrations of PCT can increase up to 1,000 ng/mL in patients with sepsis, severe sepsis or septic shock. It is detectable within 2 to 4 hours and peaks within 6 to 24 hours (as opposed to CRP which begins to rise after 12 to 24 hours and peaks at 48 hours).^[1-3] Thus, adequate antibiotic treatment leads to decreasing PCT levels.^[4] PCT has been used as helpful diagnostic tool for sepsis and septic shock.^[5] PCT measurements are used for early diagnosis and help in the early treatment.

3. Principle of the assay

The FREND PCT is a rapid quantitative "sandwich" immunoassay using fluorescent nanoparticles which measures the concentration of PCT. A 35µL drop of patient serum is placed in the FREND PCT cartridge inlet port, where the sample interacts with a proprietary mix of dry-loaded reagents. One of these reagents includes antibody-conjugated fluorescent nanoparticles, forming immune complexes with PCT in the patient sample. Capillary action moves the sample to the detection region, where capture antibodies grab the PCT-nanoparticle. The concentration of PCT is calculated by the FREND™ System when the ratio of Test/Reference fluorescence in an unknown is compared to that same ratio for standards of known concentration. The result calculated using information stored on the lot specific FREND PCT Code chip and then is displayed on the FREND™ System screen. A hard copy printout can be obtained if desired. A ratio calculated between the Reference zone and the Test zone corrects for test-to-test variations.

Total PCT concentration in a sample analyzed with the FREND PCT on the FREND™ System correlates directly with the fluorescence intensity-the higher the PCT concentration, the greater the fluorescence ratio. The FREND PCT has a measuring range determined as 0.07ng/mL to 32.00ng/mL.

4. Material provided

Q'ty	Contents	Catalogue number
25	FREND PCT cartridge(s)	FRPC 025
30	Disposable pipette tip(s)	
01	FREND PCT Code chip	
01	FREND PCT Package Insert	

5. Materials required but not provided

- The FREND™ System
- Micro-pipette or an equivalent pipette capable of delivering 35 and 70 µL
- Personal protective equipment and biohazard waste disposal containers

6. Warning and Precautions

- PCT cartridges are only to be used on the FREND™ system.
- PCT cartridges are disposable, single use devices. Do not reuse them under any circumstances.
- Allow sealed cartridges to come to room temperature for 15-30 minutes prior to use.
- Cartridges should not be frozen.
- Assure the humidity in the laboratory is in the 10-80% range when tests are run.
- Assure the room temperature remains in the range of 22-30°C when tests are run.
- Avoid cross-contamination between samples by using a new pipette tip for each new Specimen.
- Avoid high humidity, direct sunlight or heat in the area used for cartridge storage.
- Inaccurate results are possible if the sample used is contaminated in any way.
- Using specimens containing clotted fibrin could result in erroneous results.
- Over or under loading the cartridge with sample may result in inaccurate results.
- Human specimens are not used in the preparation of this product, however, since human specimens will be used for samples and other quality control products in the lab may be derived from human materials. Please use Universal Precautions when handling all specimens and controls.
- Do not use the cartridges beyond the expiration date on the pouch
- Do not use the cartridge if the pouch is damaged or the seal is broken.
- Perform testing as specified in the Package Insert and User Manual.
- Keep the cartridge sealed in the pouch until just ready for use.
- Use the cartridge immediately after opening the pouch.
- Wear disposable gloves when handling the cartridges and the samples.
- Wash hands thoroughly and often after handling reagent cartridges or samples.
- Do not ingest the silica gel package found in the cartridge pouch.
- PCT has been designed so that the high dose "hook effect" does not affect the vast majority of samples.
- Handle specimens in accordance with the OSHA Standard on Bloodborne Pathogens.

7. Storage and Stability

All unopened materials are stable until the expiration date on the label when stored at the specified temperature. Reagent stability has been demonstrated for eighteen months from the date of manufacture.

The expiration date is clearly indicated on the product box and the cartridges.

Materials	Catalogue number
Refrigerator temperature storage (2~8 °C)	
FREND PCT cartridges	FRPC 025
Room temperature storage (18~25 °C)	
Pipette tips	None

8. Specimen collection and handling

Human serum samples are suitable for use with FREND PCT cartridges. Follow instructions detailed in this package insert as well as the specimen collection tube manufacturer's instructions for specimen collection and preparation (including manufacturer's instructions for centrifugation time and speed).

For serum, a blood sample is collected aseptically without additives by venous puncture. After allowing the sample to clot for 30 minutes at room temperature, the collection tube should be centrifuged for 10 minutes at 3,000rpm. Samples may be stored at 2-8 °C for up to 6 hours prior to analysis. If the analysis is scheduled to be done at some later time, the sample should be stored frozen at -20 °C or below for future use. Repeated freeze-thaw cycles should be avoided. Prior to assay, slowly bring frozen samples to room temperature (18-25 °C) and mix gently but thoroughly before test. For optimal results, avoid grossly hemolytic, lipemic, or turbid specimens. Specimens should be free of aggregated fibrin, red blood cells, or other particulate matter.

When pipetting into the FREND PCT cartridge sample inlet, ensure that bubbles in the sample are avoided. Bubbles may restrict flow and result in an incomplete or erroneous test result.

9. Procedure

Calibration

There is no need for calibration to be performed by the end user as is generally required on other automated laboratory equipment. All calibration statistics and information have been electronically stored on the FREN D PCT Code chip included in each box of FREN D PCT cartridges. The FREN D PCT Code chip is specific for each manufactured lot of FREN D PCT cartridges.

Calibration information should always be checked by running external quality control samples to verify that the results obtained for PCT on the FREN D™ System using the FREN D PCT cartridges of a specific lot meet the laboratory criterion for acceptability.

Code chip installation

Please refer to the FREN D™ System user manual for more detailed instructions relative to the Code chip installation. Abbreviated instructions follow here:

- 1) Insert the FREN D™ System electrical cord into an appropriate outlet.
- 2) Insert the Code chip into the Code chip slot at the rear of the FREN D™ System following the arrows.
- 3) Press the **'Setup'** button on the **'Main'** screen.
- 4) Press the **'Code chip'** button on the **'Setup'** screen.
- 5) The information embedded on the FREN D PCT Code chip is automatically saved on the FREN D™ System.
- 6) When the Code chip installation is completed, press the **'OK'** button to go to the **'Setup'** screen.
- 7) Press the **'Item'** button on the **'Setup'** screen.
- 8) Check the FREN D PCT cartridge lot number and the installation date of the Code chip.
- 9) Press the **'Home'** button to go to the **'Main'** screen to begin running external quality control and patient samples.

Quality control

1. FREN D™ System QC cartridges

FREN D QC Cartridge contains multiple controls to check optic part of the system. By testing QC Cartridge, part of analytical components of the system of (1) laser power, (2) alignment, and (3) mechanical integrity are confirmed. For each day of patient testing perform QC Cartridge testing. Refer to the quality control procedures section in the User manual of FREN D™ System.

In brief, perform QC Cartridge testing for the following conditions:

- 1) Upon initial setup of the system
- 2) Each day of patient testing,
- 3) When the system has been transported or moved,
- 4) Whenever there is uncertainty about the performance of the system,
- 5) Whenever required by your laboratory's quality control requirements.

2. Internal procedural controls

FREND PCT cartridge contains built-in control features. Fluorescence signal in the Reference zone of each cartridge shows: (1) that enough volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this Reference zone signal is missing or lower than threshold, the FREND™ System consider it as an incorrect or failed test, not producing a test result but an error message. In addition, with each cartridge run, the system monitors, in part, for (1) flow of samples, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.

3. External quality control testing

Commercially available controls from a variety of manufacturers are available that contain PCT as a measured analyte. It is recommended that a minimum of two (2) levels of controls be run at least once per month or once for each new lot, whichever comes earlier. However, Controls should be run with a minimum frequency, depending on number of tests run in the laboratory. Each laboratory should establish its own criteria based on the following parameters:

- 1) Each new lot,
- 2) Each new shipment (even if from the same lot previously received),
- 3) Each new operator (an individual who has not run the tests for at least two weeks),
- 4) Monthly, as a continued check on storage conditions,
- 5) Whenever problems (storage, operator, or other) are identified,
- 6) Or other times as required by your laboratory's standard QC procedures.

Individual laboratory policy will dictate exactly which control materials and lot numbers should be run, the frequency with which controls are to be tested, criteria for acceptance of the results and required corrective action to be taken if results do not meet laboratory criteria. If any external quality control sample values are out of the acceptable range, it will be necessary to investigate the problem before reporting patient results to assure there is not

an instrument or software malfunction. Do not assay patient samples on the FREND™ System using FREND PCT if quality control results do not give expected values. Refer to your laboratory policies on how to determine acceptability of external control materials results. Each laboratory operates under a different set of regulations. Every laboratory must follow the standardized procedures acceptable to the regulatory agencies to which the laboratory is responsible.

Specimen processing

1. Preparation

Remove sufficient cartridges of FREND PCT from the refrigerator to test the number of patient samples and required external quality materials. Allow the tubes and the sealed pouches containing the cartridges to come to room temperature for 15-30 minutes prior to the start of the testing sequence.

If using refrigerated patient samples, remove those from the refrigerator and allow to them to come to room temperature prior to testing. If frozen samples will be utilized, be sure these are removed from the freezer, thawed naturally and then mixed gently but thoroughly prior to testing. Testing should not begin on these previously frozen samples until they have reached room temperature.

There are no other reagents or sample preparations necessary.

2. Assay procedure

- 1) Prepare the FREND PCT cartridge and specimen.
- 2) Pipette 35 µL of the specimen into the sample inlet on the cartridge using a suitable micro-pipette equipped with a fresh pipette tip.
- 3) Press the 'Test' button on the 'Main' screen of the FREND™ System.
- 4) The system moves to the Patient ID screen automatically.
- 5) Type the Patient ID and press the 'Enter' button to begin the test.
- 6) Insert the cartridge into the cartridge slot using the cartridge arrow as a guide.

 **Caution:** Please check the direction of the cartridge before insertion and assure the insertion is complete. It is recommended to insert the cartridge after the sample injection after 30 seconds elapsed in less than 5 minutes to obtain the optimal result of test.

- 7) When the reaction in the cartridge is completed, the FREND™ System will automatically begin the reading process.
- 8) When the measurements are completed, the cartridge will automatically be expelled and the results displayed.

 **Caution:** Do not disconnect power cord or shut off power from the FREND™ System while a cartridge is in the reading chamber. This may cause a system error.

- 9) If the FREND™ System is connected to the optional printer, press the 'Print' button and the results will be output on the printer paper.
- 10) For more detailed instructions, please refer to the FREND™ System User manual.

10. Procedural notes

If a specimen Procalcitonin concentration is found to be greater than the linearity limit of the assay of 32.00 ng/mL and a definitive result is required, the specimen should be diluted with low concentration sample that has been previously measured on the FREND PCT and found to contain < 0.07 ng/mL PCT and then re-assayed according to the 'Assay Procedure'. The recommended dilution for samples with an initial result of >32.00 ng/mL is 1:10. It is desirable to dilute the sample so that the diluted sample reads between 2.5 and 30 ng/mL. Dilutions must be made manually and the final result on the diluted sample calculated manually by multiplying the result obtained on the diluted sample by the dilution factor.

※ **Original Unknown Concentration of PCT = Concentration of diluted sample x 10**

11. Calculation of results

The FREND™ System performs all sample and reagent handling operations automatically within the cartridge once the sample has been manually loaded to the sample inlet in the cartridge and the cartridge placed into the FREND™ System. The rate of fluorescence produced by the reaction is read at various intervals during the analysis process, blank reading are subtracted after which the net rate is automatically converted to PCT concentration in ng/mL based upon information stored on the FREND PCT Code chip.

This result is then output on the screen and to the optional printer. It is also stored in memory on the FREND™ System.

Screen displayed for various concentration scenarios

Displayed result	Description
	PCT Concentration Less than 0.07 ng/mL
	PCT Concentration Not less than 0.07 ng/mL And not higher than 32.00 ng/mL
	PCT Concentration Higher than 32.00 ng/mL

12. Limitations of the procedure

- 1) When used for diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g. symptoms, results of other tests, clinical impressions, medical history, therapy, etc.)
- 2) The FRENDS™ System paired with a FRENDS PCT cartridge, is programmed to report 32.00 ng/mL as the highest concentration of PCT measurable without dilution. The lowest measurable concentration is 0.07 ng/mL - the assay limit of detection

- 3) Specimens from patients with heterophilic antibodies, such as anti-mouse (HAMA), anti-goat (HAGA), or anti-rabbit (HARA) antibodies, maybe show falsely elevated or depressed values or may result in an incomplete test. Patients routinely exposed to animals or animal serum products can be prone to these types of heterophilic interferences.
- 4) Certain medications may interfere with assay performance. All results should be interpreted with respect to the clinical picture of the patient.
- 5) Although hemolysis has an insignificant effect on the assay, hemolyzed samples may indicate mistreatment of a specimen prior to assay and results should be interpreted with caution.
- 6) Lipemia has an insignificant effect on the assay except in the case of gross lipemia where interference with the lateral flow of the sample in the cartridge may occur.
- 7) The concentration of PCT in a given sample determined with assays from different manufacturers can vary due to differences in assay methods, calibration, and antibody specificity.
- 8) Please refer to the Specimen Collection and Handling, Warnings and Precautions, Storage and Stability, and Procedural Notes sections in this insert sheet.
- 9) FREND PCT has not been validated in point-of-care settings.
- 10) FREND PCT is to be used in licensed clinical laboratories with trained technicians.

13. Performance characteristics

Precision

1. Repeatability

CLSI document EP5-A2 was utilized as a guideline for precision studies on FREND PCT. Three human serums were assayed in duplicates, twice per day over a period of 20 days. The results met with acceptance criteria in repeatability test of FREND PCT.

Sample	Conc. (ng/mL)	Repeatability		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
Low	0.50	0.074	14.6	0.019	3.8	0.003	0.6	0.077	15.1
Med	2.00	0.193	9.3	0.037	1.8	0.043	2.1	0.202	9.7
High	15.00	0.672	4.5	0.151	1.0	0.273	1.8	0.741	4.9

2. Reproducibility

The results met with acceptance criteria in reproducibility test of FREND PCT.

1) Between-lot

Sample	Conc. (ng/mL)	606909	606910	606911	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	0.50	0.48933	0.49800	0.45133	0.47956	0.08279	17.26
Med	2.00	2.09267	1.95800	1.98000	2.01022	0.27338	13.60
High	15.00	14.31533	15.22000	14.18333	14.57289	0.95113	6.53

2) Between-operator

Sample	Conc. (ng/mL)	Tester 1	Tester 2	Tester 3	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	0.50	0.46000	0.52000	0.46667	0.48222	0.07661	15.89
Med	2.00	2.14000	2.15333	2.22000	2.17111	0.26303	12.12
High	15.00	14.70667	14.36333	15.29667	14.78889	0.81415	5.51

3) Between-site

Sample	Conc. (ng/mL)	Tester 1 Lab A Instrument 1	Tester 2 Lab B Instrument 2	Tester 3 Lab C Instrument 3	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	0.50	0.54133	0.52333	0.57733	0.54733	0.09616	17.57
Med	2.00	2.02200	2.21867	2.04533	2.09533	0.29579	14.12
High	15.00	14.66267	14.05067	15.07000	14.59444	0.90973	6.23

Accuracy


The Accuracy was performed using two concentrations of PCT standard samples. The results were measured within criteria ($85\% \leq \text{Recovery} \leq 115\%$).

Comparative analysis

FREND PCT was compared to the predicate device using guidelines outlined in CLSI document EP09-A2. Samples (n=102) were measured in duplicate on each systems. Linear regression analysis demonstrated a correlation coefficient (r) of ≥ 0.95 .


1. Predicate device A

Slope	Correlation coefficient (r)
1.0014	0.9877

2. Predicate device B

Slope	Correlation coefficient (r)
2.0062	0.9849

Analytical sensitivity

The limit of blank (LoB) and limit of detection (LoD) was determined using guidelines found in CLSI document EP17-A. LoB was determined using 20 replicate measurements of five PCT negative patient samples (0.00 ng/mL). LoD was determined using 20 replicates measurements of five low patient samples.

LoB	LoD	LoQ
0.04 ng/mL	0.07 ng/mL	0.07 ng/mL

Linearity

The range of linearity in serum was established using CLSI document EP6-A. FERND PCT was found to have linearity within the reportable range of 0.07 ng/mL ~ 39.10 ng/mL. The dilution linearity study was performed by diluting a high concentration PCT specimen with a low concentration PCT specimen. The analytical measurement range of FRENED PCT is from 0.07 ng/mL to 32.00 ng/mL.

Cross-reactivity

The following substances were evaluated for potential cross-reactivity with the FRENED PCT at three concentrations. Testing was done according to the CLSI protocol EP7-A2. No significant cross-reactivity was found.

No.	Substances	Concentration
1	Albumin	4 g/dL
2	Human calcitonin	60 ng/mL
3	Human Katakalcin	10 ng/mL
4	Human a-cGRP*	10 µg/mL
5	Human b-cGRP*	10 µg/mL

**Calcitonin Gene Related Peptide*

Interference

1. Endogenous interference

The interference study was performed as recommended in the CLSI EP7-A2 protocol using three concentrations of PCT. No interference by the substances below was found.

No.	Substances	Concentration
1	Hemoglobin	500 mg/dL
2	Bilirubin	40 mg/dL
3	Triglycerides	22.5 mg/mL

2. Pharmaceutical interference

A variety of common medications were tested for interference with the FREND PCT. The testing showed there was no significant interference (<10%) from the tested drugs that would affect the interpretation of a procalcitonin result as assayed on the FREND PCT.

No.	Substances	Concentration
1	Imipenem	0.5 mg/mL
2	Cefotaxime	180 mg/dL
3	Vancomycin	3 mg/mL
4	Dopamine	26 mg/dL
5	Noradrenalin	4 µg/mL
6	Dobutamin	22.4 µg/mL
7	Heparin	16,000 U/L
8	Furosemide	4 mg/dL

14. References

- 1) Christ-crain and Muller. Biomarkers in respiratory tract infections: diagnostic guides to antibiotic prescription, prognostic markers and mediators. *Eur Respir J.* 2007, 30:556-73.
- 2) Kibe S, et al. Diagnostic and prognostic biomarkers of sepsis in critical care. *J Antimicrob Chemother.* 2011, 66 (S2). 33-40.
- 3) Schuetz P, et al. Procalcitonin for diagnosis of infection and guide to antibiotic decisions: past, present and future. *BMC Medicine.* 2011, 9: 107.
- 4) Gattas and Cook, Procalcitonin as a Diagnostic Test for Sepsis: Health Technology Assessment in the ICU. *J Crit Care.* 2003. 18 : 52-8
- 5) Becker KL, et al. Procalcitonin assay in systemic inflammation, infection, and sepsis: clinical utility and imitations. *Crit Care Med.* 2008 March, 36 (3): 941-52

Glossary of Symbols

	Caution, warning, Consult accompanying documents
	Catalogue number/Reference number
 www.nanoentek.com/eifu.php	Consult Instructions for Use An electronic instructions for use (eIFU) indicator (website address) may accompany the symbol when used to indicate an instruction to consult an eIFU.
	Lot number/Batch number
	Use by YYYY-MM-DD or YYYY-MM
	Manufacturer
	CE marking
	<i>In vitro</i> diagnostic medical device
	Temperature limitation
	Contains sufficient for <n> tests
	Do not reuse
	Do not use if package is damaged
	For prescription use only CAUTION: Federal (U.S.) law restricts this device to sale by or on order of a physician.
	US Corporation
	Patient ID
	Result
 Sample Drop	Sample Drop
	Authorized representative in the European Community

UK Representative	Authorized representative in United Kingdom
CH REP	Authorized representative in Switzerland
BRH	Authorized representative in Brazil

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NanoEntek, Inc.
MT Promedt Consulting GmbH