BUDDI™ Influenza A & B Rapid Card

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BUDDI™ Influenza A & B Rapid Card

※ Product code : 11371

※ Please read this Instructions before use.

INTENDED PURPOSE

It is In vitro diagnostic medical device to help diagnosis of influenza A virus and influenza B virus by qualitative test from nasal or throat swab specimens of the patients.

EXPLANATION OF THE TEST

Human influenza A and B viruses are major infectious pathogens of the respiratory tract that affect people of all ages, with the highest rates of mortality and severe morbidity in the elderly and the young. Patients who present with suspected influenza may benefit from treatment with an antiviral agent especially if given within the first 48 hrs of onset of illness. It is important to distinguish influenza A from influenza B in order to allow physicians a choice in selective antiviral intervention. Diagnostic tests available for influenza include rapid immunoassay, immunofluorescence assay, polymerase chain reaction(PCR), serology, and viral culture.

BUDDI™ Influenza A&B Rapid Card is a qualitative rapid test to detect influenza A or B antigens from respiratory specimens of symptomatic patients

PRINCIPLE OF THE METHOD

BUDDI™ Influenza A&B Rapid Card is based on the principle of an immunochromatography in vitro test for the qualitative determination of influenza A or influenza B. When the test device including specimens dissolved in extraction solution, the specimens moves to the conjugate pad via sample pad and resuspends anti-influenza A or anti-influenza B conjugated gold complex that is dried on the conjugate pad. The mixture moves along the membrane by capillary action and reacts with anti-influenza A or anti-influenza B that is immobilized on membrane. If influenza A or influenza B is present enough in the sample, a colored band is appeared. If there is no influenza A or influenza B in the sample or not sufficient, the area will remain colorless. As the control band(C) which is red or purple color is appeared, it means the test is working properly and the result is valid.

COMPONENTS

- ① Test device ② Specimen extraction tube ③ Filter cap ④ Specimen collection swab ⑤ Instructions for use


*Control swabs are optional.

SPECIMEN COLLECTION AND PREPARATION

01. Specimen collection

(01) Nasal swab specimens

To collect nasal swab specimens, insert the sterile swab into the nostril that presents the most secretion under visual inspection. Using gentle rotation push the swab until resistance is met at the level of the turbinate(less than one inch into the nostril). Rotate the swab a few times against the nasal wall.

Nasal swab

Throat swab

(02) Throat swab specimens

To collect throat swab specimens, rub the sterile swab on both tonsillar surfaces and the posterior pharynx.


02. Specimen transport and storage

Testing should be performed immediately after the specimens have been collected. The specimens extracted from swab may be stored at 1-8°C for up to one day. For long term storage, specimens should be kept below -20°C.

TEST PROCEDURE

- All reagents, specimens and test devices must be at room temperature (15-30°C).
- Insert the patient sample swab into extraction tube and swirl the swab 10 times.
- Remove the swab while pressing the swab head against tests tube to extract the liquid from the swab.
- Discard the used swab.
- Close the tube with filter cap.
- Take out test device from pouch and place it on a flat surface.
- Add 3 drops into the sample well(S) of the test device.
- After adding sample, reaction will be finished within 8 minutes.
- Insert test device into the BUDDI™ and result will be displayed when the reading is finished.

<Test Procedure>


READING AND INTERPRETATION OF RESULT

It is recommended to using Influenza BUDDI™ Reader.

01. NEGATIVE


If the Reader display 'Negative', It means Negative for the present of Influenza A and Influenza B antigen.

02. POSITIVE


(1) Positive for Influenza virus A

If the Reader display 'Positive' for Influenza A and 'Negative' for Influenza B, it means positive for the present of Influenza A antigen.


(2) Positive for Influenza virus B

If the Reader display 'Positive' for Influenza B and 'Negative' for Influenza A, it means positive for the present of Influenza B antigen.

(3) Positive for Influenza virus A&B

If the Reader display 'Positive' for Influenza A and 'Positive' for Influenza B, it means positive for the present of Influenza A&B antigen.

03. INVALID

If the control line ('C' band) does not present in the result window of test device, re-test is needed.

FOLLOW-UP ACTION

- BUDDI™ Influenza A&B Rapid Card** is a qualitative test kit. For knowing the quantitative concentration of Influenza, it should be performed with additional quantitative test.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

LIMITATIONS OF THE METHOD

- BUDDI™ Influenza A&B Rapid Card** is for the qualitative detection of influenza A and B viral antigens. The test performance depends on antigen load and may not correlate with cell culture performed on the same specimen. Negative test results are not intended to rule out other non-influenza viral infections.
- Sensitivity can differ with various strains of influenza due to difference in antigen expression. Specimens might contain new, non-identified strains of influenza that express varying amounts of antigen.
- As with all diagnostic tests, the results obtained with this test kit yield data that must be used only as an adjust to other information available to the physician.
- Individuals who received nasally administered influenza vaccine may have positive test results for up to three days after vaccination.
- Monoclonal antibodies may fail to detect(or with less sensitivity) influenza A viruses that have undergone minor amino acid changes in the target epitope region.

PERFORMANCE CHARACTERISTICS

01. Sensitivity and Specificity

A study was performed using 400 positive and negative specimens. Each specimen was assayed with the **BUDDI™ Influenza A&B Rapid Card** and a commercially available DFA after cell culture according to the respective instructions for use. Compared to culture and confirmed for influenza A or influenza B by DFA, the **BUDDI™ Influenza A&B Rapid Card** correctly identified 95.5 %(191/200) positive specimens and 100 %(200/200) negative specimens.

Cell Culture Test	Influenza A&B Rapid Card			Total	
		A	B		Negative
	A	97	-		5
B	-	94	4	98	
Neg.	-	-	200	200	
Total	97	94	209	400	

※ Sensitivity : Influenza A 95.10 %(97/102), Influenza B 95.92 %(94/98)

※ Specificity : more than 99.9%(200/200)

02. Analytical specificity and cross reactivity

BUDDI™ Influenza A&B Rapid Card was evaluated with 8 bacterial and 8 viral isolates. These studies had been performed to demonstrate that have been no cross reactivity with 8 bacterial panels and 8 virus panels on the **BUDDI™ Influenza A&B Rapid Card**.

03. Analytical sensitivity

Analytical sensitivity was established using a total 8 strains of influenza virus : 6 influenza A and 2 influenza B.

04. Precision

- Within run performance test was determined by one analyst with ten devices each of 3 different lots containing different concentration every influenza type. There was no variation between the test devices of each lot and between three different lots.
- Between run performance test was determined by three analysts with 3 different lots containing different concentration every influenza type. There was no variation between different analysts.

05. Analytical Interference test

The following potential interfering substances were tested and were found to have no affect on the performance of the **BUDDI™ Influenza A&B Rapid Card**.

Potential interfering substances	Concentration	Potential interfering substances	Concentration
Acetyl Salicylic Acid	20 mg/mL	Phenylephrine HCl	1 mg/mL
Acetamidophenol	10 mg/mL	Fexofenadine	500 mg/mL
Chlorpheniramine	5 mg/mL	Amantadine	500 mg/mL
Diphenhydramine HCl	5 mg/mL	Ribavirin	500 mg/mL
Guaiacol Glyceryl Ether	20 mg/mL	Pseudoephedrine HCl	20 mg/mL
Oxymetazoline HCl	0.05 mg/mL	Ibuprofen	10 mg/mL

WARNING AND PRECAUTIONS

- For in vitro diagnostic use only.
- Do not eat or smoke while handling specimens.
- Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
- Do not use test kit if the packing is damaged or the seal is broken.
- Avoid splashing or aerosol formation.
- Clean up spills thoroughly using an appropriate disinfectant.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a bio hazard container.
- Do not re-use the test device.
- Use the only swab provided (Do not use the swab not provided).
- Please refer to the user manual of BUDDI™ Reader when further information is needed to prepare the test with BUDDI™ Reader.
 - Use the device with the extraction tube in same box.

STORAGE AND SHELF-LIFE

BUDDI™ Influenza A&B Rapid Card should be stored at 2-30 °C (36-86 °F). The test device is sensitive to humidity as well as to heat. Do not use it beyond the expiration date, 24 months from manufacturing.

LITERATURE REFERENCES

- Barenfanger, J., C. Drake, T. Mueller, T. Troutt, J. O'Brien, and K. Guttman. R-Mix cells are faster, at least as sensitive and marginally more costly than conventional cell lines for the detection of respiratory viruses. *J. Clin. Virol.* 2001; 22:101-110.
- Glezen W. Emerging infections: pandemic influenza. *Epidemiol. Rev.* 1996;18:64-72.
- Hayden FG, Treanor JJ, Betts RF, Lobo M, Esinhart JD, Hussey EK. Safety and efficacy of the neuraminidase inhibitor GG167 in experimental human influenza. *J. Am. Med. Assoc.* 1996;275:295-9.

Manufactured by

NanoEnTek, Inc. | 851-14, Seohae-ro, Paltan-myeon, Hwaseong-si, Gyeonggi-do, 18531, Korea
Tel : +82-2-6220-7942 | Fax : +82-2-6220-7999 | ivdst@nanoentek.com | www.nanoentek.com

MT Promedt Consulting GmbH Altenhofstrasse 80, 66386 St. Ingbert, Germany