
FREND D-Dimer

 REF FRDD 025

- English
- German
- Italian
- Spanish
- French
- Portuguese

FREND D-Dimer

Quantitative Assay for Fibrin Degradation product

REF FRDD 025

IVD For *in vitro* diagnostic use only

1. Intended use

The FREND D-Dimer is designed for *in vitro* quantitative measurement of Fibrin degradation products (FbDP, D-Dimer) in human plasma (Sodium citrate) by fluorescence immunoassay. The FREND D-Dimer is indicated for use in conjunction with a clinical pretest probability assessment model to exclude deep vein thrombosis (DVT) and pulmonary embolism (PE) disease in outpatients suspected of DVT or PE. For professional use only.

2. Principle of the assay

D-Dimer is a final product of fibrin degradation resulting from the fibrinolysis of cross-linked fibrin in blood clots by plasmin.^[1] D-Dimer consists of two D fragments of the fibrin protein linked together through covalent bonds between γ -chains and is considered a marker of fibrinolysis of fibrin clots.^[2] Elevated levels of D-Dimer are not specific to any single disease but can indicate the presence of active clotting and fibrinolysis. It is associated with abnormal blood clotting conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE), and disseminated intravascular coagulation (DIC).^[3]

In clinical practice, a D-Dimer test is often used to rule out conditions like DVT and PE, particularly in patients with low to moderate pre-test probability.^[4, 5] A negative D-Dimer test, indicating low levels of the protein,

can exclude the presence of significant clotting, reducing the need for more invasive or expensive tests like imaging studies.^[6] However, since elevated D-Dimer levels can also be seen in conditions like infection, inflammation, and cancer, it is used in conjunction with other clinical assessments to guide diagnosis and treatment.^[7]

The FREND D-Dimer test cartridge is a one-time-use rapid "Sandwich" immunoassay utilizing fluorescent nanoparticle in microfluidic flow to capture and quantify total D-Dimer in sodium citrate plasma specimens. The drop of 35 μ L patient sample is added to the FREND D-Dimer cartridge, interacting with D-Dimer antibodies conjugated fluorescent nanoparticles. The mixture moves via Capillary action to the detection region, where fluorescent nanoparticle complexes are grabbed. The fluorescence intensities from the complexes are measured and D-Dimer concentration is calculated by the FREND™ System.

3. Material provided

Q'ty	Contents	Catalogue number
25	FREND D-Dimer cartridge(s)	FRDD 025
30	Disposable pipette tip(s)	
1	FREND D-Dimer Code chip	
1	FREND D-Dimer Package Insert	

4. Materials required but not provided

- The FREND™ System
- Micro-pipette, or any pipettor, capable of delivering 35 μ L.
- Personal protective equipment and biohazard waste disposal containers

5. Warning and Precautions

- The FREND D-Dimer cartridges are intended for in vitro diagnostic use only.
- The FREND D-Dimer cartridges are only to be used on the NanoEntek FREND™ System.
- The FREND D-Dimer cartridges are disposable, single use devices. Do not reuse them under any circumstances.
- Allow sealed cartridges to come to room temperature for 15-30 minutes prior to use.

- Cartridges should not be frozen.
- Assure the humidity in the laboratory is in the 10-80% range when tests are run.
- Avoid high humidity, direct sunlight or heat in the area used for cartridge storage.
- Avoid cross-contamination between samples by using a new pipette tip for each new specimen.
- Inaccurate results are possible if the sample used is contaminated in any way.
- Using specimens containing clotted fibrin could result in erroneous result.
- Over or under loading the cartridge with sample may result in inaccurate results.
- Do not use the cartridges beyond the expiration date on the pouch.
- Do not use the cartridge if the pouch is damaged or the seal is broken.
- Perform testing as specified in the Package Insert Sheet and User Manual.
- Keep the cartridge sealed in the pouch until ready for use.
- Use the cartridge immediately after opening the pouch.
- Handle specimens in accordance with the OSHA Standard on Boodborne Pathogens.
- Human specimens are not used in the preparation of this product, however, since human specimens will be used for samples and other quality control products in the laboratory may be derived from human materials, Use Universal Precautions when handling all specimens and controls. Wear disposable gloves when handling the cartridges and the samples.
- Wash hands thoroughly and often after handling reagent cartridges or samples.
- Do not ingest the silica gel packet found in the cartridge pouch.

6. Storage and Stability

All unopened materials are stable until the expiration date on the label when stored at the specified temperature. Reagent stability has been demonstrated for twelve months from the date of manufacture.

The expiration date is clearly indicated on the product box and the cartridges.

Materials	Catalogue number
Refrigerator temperature storage (2~8 °C)	
FREND D-Dimer cartridges	FRDD 025
Room temperature storage (18~25 °C)	
Pipette tips	None

7. Specimen collection and handling

Human plasma (Sodium citrate) samples are suitable for use with FREND D-Dimer cartridges.

Follow instructions detailed in this package insert as well as the specimen collection tube manufacturer's instructions for specimen collection and preparation (including manufacturer's instructions for centrifugation time and speed.)

For Sodium citrate plasma, a venous blood sample is collected aseptically with the designated additive. After allowing the specimen to sufficiently mix with anticoagulant at room temperature, the sample tube can be centrifuged for 10 minutes at 3,000 rpm. Sample may be stored at 2~8°C for up to 6 hours prior to analysis. If the analysis is scheduled to be done at some later time, the sample should be stored frozen at -20°C or below for future use. Repeated freeze-thaw cycles should be avoided. Prior to assay, slowly bring frozen samples to room temperature (64~77°F or 18~25°C) and mix gently but thoroughly before test.

Sample may be stored at 2~8°C for up to 6 hours prior to analysis. If the analysis is scheduled to be done at some later time, the sample should be stored frozen at -20°C or below for future use. Repeated freeze-thaw cycles should be avoided. Prior to assay, slowly bring frozen samples to room temperature (64~77°F or 18~25°C) and mix gently but thoroughly before test.

For optimal results, avoid grossly hemolytic, lipemic, or turbid specimens. Specimens should be free of aggregated fibrin, red blood cells, or other particulate matter. When pipetting into the FREND D-Dimer cartridge sample inlet, ensure that bubbles in the sample are avoided. Bubbles may restrict flow and result in an incomplete or erroneous test result.

8. Procedure

Calibration

There is no need for calibration to be performed by the end user. All calibration statistics and information have been electronically stored on the FREND D-Dimer Code chip included in each box of FREND D-Dimer cartridges. The FREND D-Dimer Code chip is specific for each lot of FREND D-Dimer cartridges. Always run external quality control samples to verify that the D-Dimer results obtained on the FRENDSM System meet the laboratory criteria for acceptability for each lot of FREND D-Dimer cartridge.

Code chip installation

Please refer to the FRENDSM System User manual for more detailed instructions relative to the Code chip installation. Abbreviated instructions are as follows:

- 1) Insert the FRENDSM System electrical cord into an appropriate outlet.
- 2) Insert the Code chip into the Code chip slot at the rear of the system following the arrows.
- 3) Press the **'Setup'** button on the **'Main'** screen.
- 4) Press the **'Code chip'** button on the **'Setup'** screen.
- 5) The information embedded on the FREND D-Dimer Code chip is automatically saved on the FRENDSM System.
- 6) When the Code chip installation is completed, press the **'OK'** button to go to the **'Setup'** screen.
- 7) Press the **'Item'** button on the **'Setup'** screen.
- 8) Check the FREND D-Dimer cartridge lot number and the installation date of the Code chip.
- 9) Press the **'Home'** button to go to the **'Main'** screen to begin running external quality control and patient samples.

Quality control

1. FRENDSM System QC cartridges

The FRENDSM QC cartridge contains multiple controls to check the optics of the system. By testing with the QC cartridge, the analytical components of the system; (1) laser power (2) alignment, and (3) mechanical integrity are confirmed.

For each day of patient testing, perform QC Cartridge testing. Refer to the Quality Control Procedures section in the User manual of FREND™ System. In brief, perform QC Cartridge testing for the following conditions:

- 1) Upon initial setup of the system
- 2) Each day of patient testing,
- 3) When the system has been transported or moved
- 4) Whenever there is uncertainty about the performance of the system,
- 5) Whenever required by your laboratory's quality control requirements.

2. Internal procedural controls

The FREND D-Dimer test cartridge contains built-in control features. Fluorescence signal in the Reference Zone of each cartridge shows: (1) that enough sample volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this Reference Zone signal is missing or lower than the threshold, the FREND™ System considers it as an incorrect or failed test, and produces an error message instead of a test result. In addition, with each cartridge run, the system monitors, in part, for (1) flow of sample, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.

3. External quality control testing

Commercially available controls from a variety of manufacturers are available that contain D-Dimer as a measure analyst. It is recommended that a minimum of two (2) levels of controls be run at least once per month or once for each new lot, whichever comes earlier. However, controls should be run according to the local requirements for each laboratory. Each laboratory should establish its own criteria based on the following parameters.

- 1) Each new lot,
- 2) Each new shipment (even if from the same lot previously received),
- 3) Each new operator (an individual who has not run the tests for at least two weeks),
- 4) Monthly, as a continued check on storage conditions,
- 5) Whenever problems (storage, operator, or other) are identified,
- 6) Or other times as required by your laboratory's standard QC procedures.

Individual laboratory policy will dictate exactly which control materials and lot numbers should be run, the frequency with which controls are to be tested, criteria for acceptance if the results and required corrective action to be taken

if results do not meet laboratory criteria. If any external quality control sample values are out of the acceptable range, it will be necessary to investigate the problem before reporting patient results to assure there is not an instrument or software malfunction. Do not assay patient samples on the FREND™ System using the FREND D-Dimer if quality control results do not fall within the acceptable ranges. Each laboratory operates under a different set of regulations. Every laboratory must follow the standardized procedures acceptable to the regulatory agencies to whom the laboratory is responsible.

Specimen processing

1. Preparation

Remove sufficient FREND D-Dimer cartridges from the refrigerator to test the number of patient samples. Allow the tubes and the sealed pouches containing the cartridges to come to room temperature for 15-30 minutes prior to the start of the testing sequence.

If using refrigerated patient samples, remove those from the refrigerator and allow them to come to room temperature (64~77°F or 18~25°C) prior to testing. If frozen samples will be utilized, be sure these are removed from the freezer, thawed naturally and then mixed gently but thoroughly prior to testing. Testing should not begin on these previously frozen samples until they have reached room temperature.

There are no other reagents or sample preparations necessary.

2. Assay procedure

- 1) Prepare the FREND D-Dimer cartridges and specimen at room temperature.
- 2) Record the Sample ID on the cartridge in the designated area.
- 3) Drop the sample (35 µL) into the sample inlet on the cartridge using a calibrated micro-pipette with a fresh pipette tip.
- 4) Press the **'Test'** button on the **'Main'** screen of the FREND™ System.
- 5) The screen of FREND™ System moves to the Patient ID screen automatically. If Patient screen is displayed, Press **'Patient'** button to proceed to the Patient ID screen.
- 6) Type the Patient ID and press the **'Enter'** button to begin the test.
- 7) Insert the cartridge into the cartridge slot using the cartridge arrow as a guide.

⚠ Caution: Check the direction of the cartridge before insertion and assure the insertion is completed. It is recommended to insert the cartridge after the sample loading after 30 seconds elapsed in less than 5 minutes to obtain the optimal result of test.

- 8) When the reaction in the cartridge is completed, the FREND™ System will automatically begin the reading process.
- 9) When the measurements are completed, the cartridge will automatically be expelled and the results displayed.
⚠ Caution: *Do not disconnect power cord or shut off power on the FREND™ System while a cartridge is in the reading chamber. This may cause a system error.*
- 10) If the FREND™ System is connected to the optional printer, press the 'Print' button and the results will be output on the printer paper.
- 11) For more detailed instructions, please refer to the FREND™ System User Manual.

9. Specimen Dilution Procedures

Samples cannot be diluted for D-Dimer determinations.

Samples which read ">10,000ng/mL (FEU)" should be reported as such.

10. Calculation of results

The FREND™ System performs all sample and reagent handling operations automatically within the cartridge once the sample has been loaded to the sample inlet in the cartridge and the cartridge placed into the FREND™ System. The rate of fluorescence produced by the reaction is read at various intervals during the analysis process, blank reading is subtracted after which the net rate is automatically converted to D-Dimer concentration in ng/mL (FEU) based upon information stored on the FREND D-Dimer Code chip. This result is then displayed on the screen and can be sent to the optional printer. It is also stored in memory on the FREND™ System.

Screen displayed for various concentration scenarios

Displayed result	Description
 <p>The screenshot shows a green interface with a white box containing the following text: Order Time: 2006-07-04 08:23:11 Patient ID: Order ID: Lab ID: HANSEN1010 Item: D-Dimer Test Result: = 10.00 ng/mL Normal Range</p>	<p>D-Dimer Concentration Less than 50.00 ng/mL (FEU)</p>
 <p>The screenshot shows a green interface with a white box containing the following text: Order Time: 2006-07-04 08:23:11 Patient ID: Order ID: Lab ID: HANSEN1010 Item: D-Dimer Test Result: = 200.00 ng/mL Normal Range</p>	<p>D-Dimer Concentration Not less than 50.00 ng/mL (FEU) and not higher than 10,000.00 ng/mL (FEU)</p>
 <p>The screenshot shows a green interface with a white box containing the following text: Order Time: 2006-07-04 08:23:11 Patient ID: Order ID: Lab ID: HANSEN1010 Item: D-Dimer Test Result: > 10000.00 ng/mL Normal Range</p>	<p>D-Dimer Concentration Higher than 10,000.00 ng/mL (FEU)</p>

11. Limitations of the procedure

When used for diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g. symptoms, results of other tests, clinical impressions, medical history, therapy, etc.).

The FRENTM System paired with a FREN D-Dimer cartridge, is programmed to report 10,000 ng/ml (FEU) as the highest concentration of D-Dimer measurable without dilution. The lowest measurable concentration is 50 ng/ml (FEU) - the assay limit of quantitation. If the result is below the lowest reportable range, it should be reported as such < 50 ng/ml (FEU).

12. Performance characteristics

Within-run (Repeatability)

Clinical and Laboratory Standards Institute (CLSI) guideline EP5-A2 was utilized as a guideline for precision studies on FREN D-Dimer. It was tested twice a day repeatedly for twenty days in one lot in each sample. The results met with acceptance criteria in repeatability test of FREN D-Dimer.

Standard sample	Conc. (ng/mL)	Repeatability		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
Low	500.00	41.923	8.4	4.861	1.0	6.532	1.3	42.707	8.6
Med	2000.00	145.388	7.1	70.590	3.4	18.011	0.9	162.619	7.9
High	4000.00	313.297	7.9	155.689	3.9	82.391	2.1	359.420	9.1

Between-run (Reproducibility)

The results met with acceptance criteria in reproducibility test of FREN D-Dimer.

1) Between-lot

Standard Sample	Conc. (ng/mL)	177901	177902	177903	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	500.00	505.04	519.70	508.93	511.22	40.648	8.0
Med	2000.00	2064.91	2034.03	1971.70	2023.55	170.811	8.4
High	4000.00	4095.72	3938.44	3996.27	4010.14	315.437	7.9

2) Between-operator

Standard Sample	Conc. (ng/mL)	Tester 1	Tester 2	Tester 3	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	500.00	481.84	472.66	551.69	502.06	42.154	8.4
Med	2000.00	2073.64	2036.61	1874.11	1994.79	242.959	12.2
High	4000.00	4179.33	4091.91	3832.99	4034.74	395.653	9.8

3) Between-site

Standard Sample	Conc. (ng/mL)	Site A	Site B	Site C	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	500.00	489.49	505.38	515.99	503.62	40.68	8.1
Med	2000.00	1969.47	2050.00	1980.09	1999.85	164.076	8.2
High	4000.00	3989.37	4071.35	3984.20	4014.97	402.768	10.0


Accuracy

The Accuracy was performed using three concentrations of D-Dimer standard samples. The results were measured within criteria (85% ≤ Recovery ≤ 115%).

Linearity

The range of linearity in serum was established using CLSI document EP6-A. FERND D-Dimer was found to have linearity within the reportable range of 50.00 ng/mL ~ 10,000.00 ng/mL. The dilution linearity study was performed by diluting a high concentration D-Dimer specimen with a D-Dimer depleted serum. Results are summarized in the graph below, but may vary in individual laboratories.

FREND D-Dimer Linearity


The analytical measurement range of FREND D-Dimer is from 50.00 ng/mL to 10,000.00 ng/mL.

Method comparison

FREND D-Dimer was compared to the predicate device using guidelines outlined in CLSI document EP09-A2-IR. Samples (n=194) were measured in duplicate on both systems. Linear regression analysis demonstrated a correlation coefficient (r) =0.9764.

Slope	Correlation coefficient (R)
0.9512	0.9764

Sensitivity

The limit of blank (LoB) and limit of detection (LoD) was determined using guidelines found in CLSI document EP17-A. LoB was determined using 20 replicate measurements of five blank patient samples (<45.00 ng/mL). LoD was determined using 20 replicates measurements of five low patient samples.

LoB	LoD	LoQ
38.20 ng/mL	48.67 ng/mL	48.67 ng/mL

Interference

The interference study was performed as recommended in the CLSI EP7-A2 protocol using three concentrations of D-Dimer. No interference by the substances below was found.

No.	Substance	Concentration
1	Hemoglobin	500 mg/dL
2	Triglyceride	3 g/dL
3	Bilirubin	20 mg/dL
4	Rheumatoid factor	400 IU/mL
5	Human albumin	60 mg/mL
6	Acetaminophen	20 mg/dL
7	Acetylsalicylic acid	65 mg/dL
8	Allopurinol	4 mg/dL
9	Amikacin sulfate	10.4 mg/dL
10	Ampicillin-Na	5 mg/dL
11	Ascorbic acid	6 mg/dL
12	Atenolol	1 mg/dL
13	Caffeine	6 mg/dL
14	Captopril	0.5 mg/dL
15	Carbamazepine	3 mg/dL
16	Chloramphenicol	5 mg/dL
17	Chlorpromazine hydrochloride	0.2 mg/dL
18	Cimetidine	2 mg/dL
19	Cinnarizine	3 mg/dL
20	Creatinine	30 mg/dL
21	Cyclosporine	0.4 mg/dL
22	Dextran 75	2500 mg/dL
23	Digoxin	0.0006 mg/dL
24	D-L methyl dopa hydrochloride	1.8 mg/dL
25	Dopamine hydrochloride	0.1 mg/dL
26	Erythromycin	6 mg/dL
27	Ethanol	400 mg/dL
28	Ethosuximide	25 mg/dL
29	Furosemide	6 mg/dL
30	Gentamicin sulfate	1 mg/dL
31	Heparin Lithium salt	300 U/dL

32	Heparin sodium salt	300 U/dL
33	Ibuprofen	50 mg/dL
34	Lidocaine	1.2 mg/dL
35	Lithium chloride	14 mg/dL
36	Levothyroxine	0.06 mg/dL
37	Nicotine	0.1 mg/dL
38	Nifedipine	0.04 mg/dL
39	Penicillin G sodium salt	2500 U/dL
40	Phenytoine	5 mg/dL
41	Primidine	4 mg/dL
42	Propranolol hydrochloride	0.2 mg/dL
43	Theophylline	4 mg/dL
44	Urea	500 mg/dL
45	Uric acid	24 mg/dL
46	Valproic acid sodium salt	60 mg/dL
47	Verapamil hydrochloride	0.2 mg/dL
48	Warfarin	1 mg/dL

Cross reactivity





The following substances were evaluated for potential cross-reactivity with the FREN D-Dimer at three concentrations. Testing was done according to the CLSI protocol EP7-A2. There was no cross-reactivity in the test of the FREN D-Dimer with fibrinogen, FDP-X, FDP-D, and FDP-E.

No.	Substance	Concentration
1	Fibrinogen	1 mg/dL
2	FDP-X	10 µg/dL
3	FDP-D	1 µg/dL
4	FDP-E	10 µg/dL

13. References


- 1) Adam SS, Key NS, Greenberg CS. D-Dimer antigen: current concepts and future prospects. *Blood*. 2009 Mar 26;113(13):2878-87.
- 2) Hoeprich PD Jr, Doolittle RF. Dimeric half-molecules of human fibrinogen are joined through disulfide bonds in an antiparallel orientation. *Biochemistry*. 1983 Apr 26;22(9):2049-55
- 3) Wells PS. Integrated strategies for the diagnosis of venous thromboembolism. *J Thromb Haemost*. 2007 Jul;5 Suppl 1:41-50.
- 4) van der Hulle T, den Exter PL, Erkens PG, van Es J, Mos IC, ten Cate H, Kamphuisen PW, Hovens MM, Büller HR, Klok FA, Huisman MV. Variable D-Dimer thresholds for diagnosis of clinically suspected acute pulmonary embolism. *J Thromb Haemost*. 2013 Nov;11(11):1986-92.
- 5) van Es N, van der Hulle T, van Es J, den Exter PL, Douma RA, Goekoop RJ, Mos IC, Galipienzo J, Kamphuisen PW, Huisman MV, Klok FA, Büller HR, Bossuyt PM. Wells Rule and D-Dimer Testing to Rule Out Pulmonary Embolism: A Systematic Review and Individual-Patient Data Meta-analysis. *Ann Intern Med*. 2016 Aug 16;165(4):253-61.
- 6) Kline JA, Hogg MM, Courtney DM, Miller CD, Jones AE, Smithline HA. D-Dimer threshold increase with pretest probability unlikely for pulmonary embolism to decrease unnecessary computerized tomographic pulmonary angiography. *J Thromb Haemost*. 2012 Apr;10(4):572-81.
- 7) Weitz JI, Fredenburgh JC, Eikelboom JW. A Test in Context: D-Dimer. *J Am Coll Cardiol*. 2017 Nov 7;70(19):2411-2420.

Glossary of Symbols

	Caution, warning, Consult accompanying documents
	Catalogue number/Reference number
 www.nanoentek.com/eifu.php	Consult Instructions for Use An electronic instructions for use (eIFU) indicator (website address) may accompany the symbol when used to indicate an instruction to consult an eIFU.
	Lot number/Batch number
	Use by YYYY-MM-DD or YYYY-MM
	Manufacturer
	CE marking
	<i>In vitro</i> diagnostic medical device
	Temperature limitation
	Contains sufficient for <n> tests
	Do not reuse
	Do not use if package is damaged
	For prescription use only CAUTION: Federal (U.S.) law restricts this device to sale by or on order of a physician.
	US Corporation
	European Corporation
	Patient ID
	Result
	Sample Drop

EC REP	Authorized representative in the European Community
UK Representative	Authorized representative in United Kingdom
CH REP	Authorized representative in Switzerland
BRH	Authorized representative in Brazil

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NanoEntek, Inc.
NanoEntek Europe | med-tech supplies GmbH
MT Promedt Consulting GmbH