
FREND hsCRP

 FRhCR 025

- English
- German
- Italian
- Spanish
- French
- Portuguese

FREND hsCRP

Quantitative assay for high sensitivity C-reactive Protein

REF FR_hCR 025

IVD For *in vitro* diagnostic use only

1. Intended use

FREND hsCRP is a fluorescence immunoassay (FLA) for use with the FREND™ System. It is designed for the quantitative detection of the human C-reactive protein in human serum which may help in the evaluation of infection, tissue injury, inflammatory disorders and associated diseases. It may also aid in the identification of individuals at risk for future cardiovascular disease. In conjunction with traditional evaluation method FREND hsCRP may be useful as a prognostic marker for stable coronary disease or acute coronary syndromes. ^[1,2] For professional use only.

2. Principle of the assay

FREND hsCRP is a single-use quantitative fluorescence immunoassay (FIA) kit that can measure human CRP concentration present in human serum. When 35 µL of human serum is transferred into the cartridge inlet, the CRPs presented in human serum are captured by antibody specific to human CRP and detected by anti-CRP antibody conjugated to fluorescent microparticles.

3. Material provided

Q'ty	Contents	Catalogue number
25	FREND hsCRP cartridge(s)	FRhCR 025
30	Disposable pipette tip(s)	
1	FREND hsCRP Code chip	
1	FREND hsCRP Package Insert	

4. Materials required but not provided

- Micro-pipette or an equivalent pipette capable of delivering 35 and 70 μ L
- FREND™ System
- Personal protective equipment and biohazard waste disposal containers

5. Warning and Precautions

- FREND hsCRP cartridges are intended for in vitro diagnostic use only.
- FREND hsCRP cartridges are only to be used on the FREND™ System.
- FREND hsCRP is disposable, single use device. Do not reuse them under any circumstance.
- Allow sealed cartridges to come to room temperature for 15-30 minutes prior to use.
- Cartridges should not be frozen.
- The humidity in the laboratory must be 10-80% range when running tests.
- Avoid cross-contamination between samples by using a new pipette tip for each new specimen.
- Avoid high humidity, direct sunlight or heat when storing cartridges.
- Testing of contaminated samples may cause erroneous results.
- Inadequate pipetting may deliver under or over volumes of specimen into the cartridge which may affect test result.
- Discard and do not use any damaged or dropped cartridges.
- Do not use the cartridges beyond the expiration date on the pouch.
- Do not use the cartridge if the pouch is damaged or the seal is broken.
- Perform testing as specified in the package insert and user manual.
- Inadequate or inappropriate sample collection, storage, and transport may affect false test results.

- Keep the cartridge sealed in the pouch until ready to use.
- Use the cartridge immediately after opening the pouch.
- For professional use only.
- Use Universal Precautions when handling all specimens and controls.
- Wear disposable gloves when handling the cartridges and the samples.
- Wash hands thoroughly and often after handling reagent cartridges or samples.
- Do not ingest the silica gel packet found in the cartridge pouch.
- Do not bend cartridges.

6. Storage and Stability

All unopened materials are stable until the expiration date on the label when stored at the specified temperature. Reagent stability has been demonstrated for twelve months from the date of manufacture.

The expiration date is clearly indicated on the product outer box and the cartridges.

Materials	Catalogue number
Refrigerator temperature storage (2~8 °C)	
FREND hsCRP cartridges	FRhCR 025
Room temperature storage (18~25 °C)	
Pipette tips	None

7. Procedure

Calibration

There is no need for calibration to be performed by the end user as is generally required on other automated laboratory equipment. All calibration statistics and information have been electronically stored on the FREND hsCRP Code chip included in each box of FREND hsCRP cartridges. The FREND hsCRP Code chip is specific for each manufactured lot of FREND hsCRP cartridges.

Calibration information should always be checked by running external quality control samples to verify that the results obtained for hsCRP on the FREND™ System using the FREND hsCRP cartridges of a specific lot meet the laboratory criterion for acceptability.

Code chip installation

A lot-specific Code chip is supplied with each kit of FRENDS[™] hsCRP. When using a new lot of reagent, the Code chip of the same lot must be installed in the FRENDS[™] System. Please refer to the FRENDS[™] System User manual for more detailed instructions relative to the Code chip installation. Abbreviated instructions are as follows:

- 1) Insert the FRENDS[™] System electrical cord into an appropriate outlet.
- 2) Insert the Code chip into the Code chip slot at the rear of the FRENDS[™] System following the arrows.
- 3) Press the **'Setup'** button on the **'Main'** screen.
- 4) Press the **'Code chip'** button on the **'Setup'** screen.
- 5) The information embedded on the FRENDS[™] hsCRP Code chip is automatically saved on the FRENDS[™] System.
- 6) When the Code chip installation is completed, press the **'OK'** button to go to the **'Setup'** screen.
- 7) Press the **'Item'** button on the **'Setup'** screen.
- 8) Check the FRENDS[™] hsCRP cartridge lot number and the installation date of the Code chip.
- 9) Press the **'Home'** button to go to the **'Main'** screen to begin running external quality control and patient samples.

⚠ Caution: *If a new LOT product, remove the old code chip and apply the new code chip included in the box of the new LOT product to the FRENDS[™] System.*

Quality control

• FRENDS[™] System QC cartridges

FRENDS[™] QC Cartridge contains multiple controls to check optic part of the system. By testing QC Cartridge, part of analytical components of the system of (1) laser power, (2) alignment, and (3) mechanical integrity are confirmed.

For each day of patient testing perform QC Cartridge testing. Refer to the quality control procedures section in the User manual of FRENDS[™] System. In brief, perform QC Cartridge testing for the following conditions:

- 1) Upon initial setup of the system
- 2) Each day of patient testing,
- 3) When the system has been transported or moved,
- 4) Whenever there is uncertainty about the performance of the system,
- 5) Whenever required by your laboratory's quality control requirements.

- **Internal procedural controls**

FREND hsCRP test cartridge contains built-in control feature. Fluorescence signal in the reference zone of each cartridge shows: (1) that enough volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this reference zone signal is missing or lower than threshold, the FREND™ System consider it as an incorrect or failed test, not producing a test result but an error message. In addition, with each cartridge run, the system monitors, in part, for (1) flow of sample, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.

- **External quality control testing**

Commercially available controls from a variety of manufacturers are available that contain CRP as a measured analyte. It is recommended that a minimum of two (2) levels of controls be run once per day on days when assaying patient samples on the FREND hsCRP test.

Specimen processing

- **Preparation**

Serum is required for the assay. For serum, allow the whole blood in the collecting tube to be clotted for 30 minutes at room temperature. Centrifuge the sample for 10 minutes at 3,000 rpm within 2 hours of collection and immediately separate the serum from the packed cells.

It is recommended to use samples immediately. However, if testing is not done immediately, samples may be stored at 2-8 °C for up to 30 days prior to testing. If testing is not performed within 30 days, store at -20 °C or below. Samples can be stored frozen for up to 12 months prior to testing.

- **Assay procedure**

- 1) Allow the sample and cartridge to come to room temperature for 15~30 minutes prior to the start of the testing sequence.
- 2) Open the sealed cartridge and fill in the patient information on the cartridges.
- 3) Take 35µL of the specimen into the sample inlet on the cartridge using a pipette.
- 4) Press the 'Test' button on the 'Main' screen of the FREND™ System.
- 5) The FREND™ System moves to the patient ID screen automatically.
- 6) Type the Patient ID and press the 'Enter' button to begin the test.

- 7) Insert the cartridge into the cartridge slot using the cartridge arrows as a guide.
⚠ Caution: Please check the direction of the cartridge before insertion and assure the insertion is complete.
⚠ Caution: Please insert the cartridge into the FREND™ System after loading the sample into the cartridge.
- 8) When the reaction in the cartridge is complete, the FREND™ System will automatically begin the reading.
- 9) After the reaction ends, cartridge will be ejected and the result will be displayed and saved.
⚠ Caution: Do not remove power from the FREND™ System while a cartridge is in the reading chamber. This may cause a system error.
- 10) If the FREND™ System is connected to the printer (optional), press 'Print' button and the result will be printed out.
- 11) For more detailed instructions, please refer to the 'FREND™ System User manual'.

8. Calculation of results

The FREND™ System performs all sample and reagent handling operations automatically within the cartridge once the sample has been manually loaded to the sample inlet in the cartridge and the cartridge placed into the FREND™ System. The rate of fluorescence produced by the reaction is read at various intervals during the analysis process, blank reading are subtracted after which the net rate is automatically converted to hsCRP concentration in mg/L based upon information stored on the FREND hsCRP Code chip.

This result is then output on the screen and to the optional printer. It is also stored in memory on the FREND™ System.

Screen displayed for various concentration scenarios

Displayed result	Description
	<p>CRP Concentration Less than 0.20 mg/L</p>
	<p>CRP Concentration Higher than 0.20 mg/L, And lower than 100.00 mg/L</p>
	<p>CRP Concentration Higher than 100.00 mg/L</p>

9. Performance characteristics

- Within-run (Repeatability)**

Precision was determined as described in the CLSI guideline EP5-A2. The results met with acceptance criteria in repeatability test of the FREND hsCRP.

Sample ID	Conc. (mg/L)	Repeatability		Between-run		Between-day		Within-laboratory	
		SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)
Low	1.09	0.039	3.6	0.023	2.2	0.018	1.7	0.049	4.6
Mid1	3.29	0.210	6.5	0.212	6.6	0.113	3.5	0.319	9.9
Mid2	12.87	1.189	9.4	0.224	1.8	0.059	0.5	1.211	9.6
High	43.50	3.527	8.1	1.667	3.8	1.462	3.4	4.166	9.6

- Between-run (Reproducibility)**

The results met with acceptance criteria in reproducibility test of FREND hsCRP.

1) Between-lot

Sample ID	Conc. (mg/L)	168901	168902	168903	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	1.09	1.09	1.06	1.08	1.08	0.048	4.43
Mid1	3.29	3.23	3.21	3.29	3.25	0.284	8.74
Mid2	12.87	12.48	12.70	11.83	12.34	1.186	9.61
High	43.50	46.50	45.27	44.58	45.45	3.544	7.80

2) Between-operator

Sample ID	Conc. (mg/L)	Tester 1	Tester 2	Tester 3	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	1.09	1.10	1.09	1.09	1.09	0.017	1.55
Mid1	3.29	3.09	3.04	3.13	3.09	0.178	5.75
Mid2	12.87	12.75	11.97	12.54	12.42	0.635	5.12
High	43.50	44.26	41.75	44.22	43.41	3.293	7.59

3) Between-site

Sample ID	Conc. (mg/L)	Site A	Site B	Site C	Total Mean	SD	CV(%)
		Mean	Mean	Mean			
Low	1.09	1.07	1.08	1.10	1.08	0.051	4.74
Mid1	3.29	3.12	3.10	3.20	3.14	0.247	7.88
Mid2	12.87	12.13	11.49	11.71	11.78	0.969	8.23
High	43.50	41.90	41.15	40.93	41.32	3.706	8.97

Comparative analysis

Human serums with CRP range from 0.2 to 10.00 mg/L were tested both in the FRENDS hsCRP and the predicate device A for hsCRP and compared for the correlation. The correlation test was validated in accordance with the CLSI guideline EP9-A3.

Predicate device	Range	N	Slope	Correlation coefficient (r)
A	0.20~10.00 mg/L	123	0.9454	0.9654

Human serums with CRP range from 2.50 to 100.00 mg/L were tested both in the FRENDS hsCRP and the predicate device B for conventional CRP range and compared for the correlation. The correlation test validated in accordance with the CLSI guideline EP9-A3.

Predicate device	Range	N	Slope	Correlation coefficient (r)
B	2.50~100.00 mg/L	108	1.0057	0.9548

Analytical sensitivity

Limit of blank (LoB), Limit of detection (LoD) was determined using guidelines found in CLSI EP17-A.

LoB	LoD	LoQ
0.10 mg/L	0.18 mg/L	0.18 mg/L

Cross reactivity

The cross-reactivity study was performed as recommended in the CLSI guideline EP7-A2. It is confirmed that there is no cross-reactivity found on the following substances at its concentration.

No.	Substance	Concentration
1	Myoglobin	10 µg/mL
2	Troponin I	1 µg/mL
3	CK-MB	1 µg/mL
4	D-dimer	1,000 µg/L
5	Soluble ST2	1,000 ng/mL
6	Prostate specific antigen	100 ng/mL
7	Canine CRP	20 mg/dL




Interference

The interference study was performed as recommended in the CLSI EP7-A2 protocol.

No interference by the substances below was found.


No.	Substances	Concentration
1	Hemoglobin	2 g/L
2	Billirubin (unconjugated)	20 mg/dL
3	Billirubin (conjugated)	20 mg/dL
4	Glucose	1,000 mg/dL
5	Ascorbic acid	30 µg/mL
6	Human serum albumin	1.5 g/dL
7	Triglycerides	8.19 mg/mL
8	Rheumatoid factor	100 IU/mL
9	Cholesterol (total)	251 mg/dL

Glossary of Symbols

	Caution, warning, Consult accompanying documents
	Catalogue number/Reference number
 www.nanoentek.com/eifu.php	Consult Instructions for Use An electronic instructions for use (eIFU) indicator (website address) may accompany the symbol when used to indicate an instruction to consult an eIFU.
	Lot number/Batch number
	Use by YYYY-MM-DD or YYYY-MM
	Manufacturer
	CE marking
	<i>In vitro</i> diagnostic medical device
	Temperature limitation
	Contains sufficient for <n> tests
	Do not reuse
	Do not use if package is damaged
	For prescription use only CAUTION: Federal (U.S.) law restricts this device to sale by or on order of a physician.
	US Corporation
	Patient ID
	Result
	Sample Drop
	Authorized representative in the European Community

UK Representative	Authorized representative in United Kingdom
CH REP	Authorized representative in Switzerland
BRH	Authorized representative in Brazil

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NanoEntek, Inc.
MT Promedt Consulting GmbH