


CE

# C-Slide

*Cell Counting  
Chamber Slides*



# C-Slide

## Cell Counting Chamber Slide For Used With Automated Cell Counters

- No need for coverslips
- Disposable, eliminates the need for washing
- Reduce exposure to infectious samples and risk of contamination
- Guaranteed reproducibility and reliability
- Made of sturdy quartz grade optical plastic



Scan me



## Specifications

Sample volume : 10  $\mu$ L

Staining method : Trypan blue stain

Dimension : 75 (L) x 25 (W) x 1.8 (H) mm

## Ordering information

Catalog Number	Product Name
CS-050	C-Slide, Cell counting chamber slides 50 slides (100 tests)



website [www.nanoentek.com](http://www.nanoentek.com)  
e-mail [sales@nanoentek.com](mailto:sales@nanoentek.com)

### NanoEntek, Inc.

#### Head Office

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea  
Tel +82-2-6220-7940 / Fax +82-2-6220-7999

### NanoEntek America, Inc.

220 Bear Hill Road, Suite 102, Waltham, MA 02451, USA  
Tel +1-781-472-2558 / Fax +1-781-790-5649

### NanoEntek Europe | mts med-tech supplies GmbH

Lochhamerstr. 4a, 82152 Martinsried, Germany  
Tel +49-89-21-55-38-43 / Fax +49-89-99-95-46-60

NESCT-CS-001E (V.0.3)  
FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.